



TRUWAX

Sterile Bone Wax
Barrier to control bleeding from bone injury

Order Code :

- ◆ S 810

Description :

- ◆ Bone wax **STERILE R**

Units per box :

- ◆ 12pcs of 2.5gms each

Composition :

- ◆ Truwax is a sterile mixture of Beeswax & Paraffin wax & Isopropyl Palmitate

Indications :

- ◆ Truwax is used for the mechanical hemostasis in bone during injury / surgery.
- ◆ To patch up holes made in the skull

Most used in :

- ◆ Open heart and thoracic surgery.
- ◆ Neuro surgery
- ◆ Orthopaedic surgery

Contraindications :

- ◆ In infected areas, Where osseous regeneration and / or fusion is desired

Advantages :

- ◆ Soft, easy to form and to apply
- ◆ Easy to handle due to peel-pack with non stick cover
- ◆ Provided sterile for single use only
- ◆ Biochemically inert
- ◆ Minimally absorbable

Manufactured by :



SUTURES INDIA PVT LTD

472 D, 13th Cross, 4th Phase, Peenya Industrial Area,
Bangalore - 560 058, INDIA. Tel +91 80 4186 8000 (30 Lines)
Fax +91 80 4117 1056 Email sales@suturesin.com
www.suturesin.com