

TRULENE

Polypropylene Mesh



Non Absorbable Surgical Mesh
For Hernia Repair

Also available in ready to use Hernia Kit*



SUTURES INDIA
protecting lives

High tensile strength and burst resistance - provides permanent support.

TRULENE MESH

New design for elastic, flexible, strong, transparent Mesh.

Ideal porosity for high visibility, colonization and integration



Bio compatible, open knit construction allows in-growth of host tissues and encourages natural repair.

Maximum support and elasticity.

- No shrinkage provides long term material stability - important to counter muscular contraction.
- Thinner mesh for easy manipulation and placement for Laparoscopic surgery. No - memory suitable for any anatomic placement.

TRULENE MESH SIZES

Code	Dimensions	Packing
TM 611	6 cm x 11 cm	6 sheets / box
TM 715	7.6 cm x 15 cm	6 sheets / box
TM 151	15 cm x 15 cm	3 sheets / box
TM 303	30 cm x 30 cm	1 sheet / box



HERNIA KIT



HERNIA KIT

Code No	Mesh	Sutures
TK 611	TM 611	SN 841, SN 842, SN 2407
TK 715	TM 715	SN 841, SN 842, SN 2407

KIT CONTENTS

Code	Material	Size	Length	Needle
TM 611	Polypropylene Mesh	6cm x 11cm	—	—
TM 715	Polypropylene Mesh	7.6cm x 15cm	—	—
SN 841	Polypropylene Suture	2 - 0	70cm	½ circle round bodied 30mm
SN 842	Polypropylene Suture	0	70cm	½ circle round bodied 30mm
SN 2407	Polyglycolic Acid Suture (Undyed)	3 - 0	70cm	¾ circle reverse cutting 26mm



SUTURES INDIA PRIVATE LIMITED

472-D, 13th Cross, 4th Phase
Peenya Industrial Area, Bangalore-560 058.
Tel : 91-80-5127 2102 / 03 / 04 / 05
Fax : 91-08-5117 1056, Email : sales@suturesin.com

CE
0434

www.suturesin.com