

## Pre-Cut Sutures

LENGTH	SIZE			
	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)
150 cm	S 2612	S 2613	S 2614	S 2615
180 cm		S 2613 MDT	S 2614 MDT	S 2615 MDT

## Technical Details

- **Suture**  
TRUGLYDE Polyglycolic Acid Surgical Suture
- **Type**  
Braided & Coated
- **Material Color**  
Violet, Undyed
- **Coating Material**  
Calcium Stearate & Polycaprolactone
- **Time of Absorption**  
Absorbed by Hydrolysis, Absorption within 60-90 days
- **Range**  
USP Sizes available - 5-0 to 2.  
Standard sutures - 150cm & 180cm  
With needle - Suture lengths from 20cm to 150cm & Needle length from 12mm to 50mm
- **Packing Color**  
Violet
- **Shelf life**  
5 Years

- **Sterilization**  
Ethylene Oxide
- **Packing**  
Primary packaging  
Aluminium foil  
Secondary packaging  
Medical grade paper & Polyester laminated poly film.  
Pack of 12pcs  
400gsm Super Chromo board.  
Master shipper  
100 boxes of 12s (5 ply of 20 Kgs/cm<sup>2</sup> burst strength - 54 x 33 x 33 cm) weighing approximately 10

## Features

- Minimal tissue reaction
- Coating ensures smooth tissue passage
- Least memory and outstanding pliability




























**SUTURES INDIA PRIVATE LIMITED**

472 D, 13th Cross, 4th Phase, Peenya Industrial Area, Bangalore-560 058. INDIA  
Tel : +91 80 4186 8000 (30 Lines) ; Fax : +91 80 4117 1056, Email : sales@suturesin.com  
[www.suturesin.com](http://www.suturesin.com)

# TRUGLYDE

























## Braided & Coated Polyglycolic Acid Suture Atraumatic Needled Sutures

NEEDLE				THREAD					
	Profile	Length	Description	Length	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)
	3/8 Circle	12 mm	 Reverse Cutting	45 cm	SN 2493	SN 2494	SN 2495		
	3/8 Circle	16 mm	 Cutting	45 cm	SN 2442		SN 2471		
	1/2 Circle	16 mm	 Round Bodied	45 70 cm	SN 2303		SN 2305		
	1/2 Circle	17 mm	 Taper Cutting	45 cm			SN 2515		
	1/2 Circle	20 mm	 Round Bodied	70 cm		SN 2304	SN 2437		
	3/8 Circle	22 mm	 Cutting	45 cm			SN 2401		
	1/2 Circle	22 mm	 Cutting	90 cm			SN 2472		
	1/2 Circle	30 mm	 Round Bodied	70 90 cm			SN 2341 SN 2317	SN 2342 SN 2338	
	Profile	Length	Description	Length	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	
	1/2 Circle Double armed	22 mm	 Round Bodied	70 cm	SN 2437 DA				
	1/2 Circle	22 mm	 Cutting	45 cm	SN 2735 A Fast(Undyed)				
	1/2 Circle	23 mm	 Reverse Cutting	35 70 cm			SN 2825 A SN 2826	SN 2826 A	
	3/8 Circle	26 mm	 Reverse Cutting	70 45 cm	SN 2407 (Undyed) SN 2732 A Fast (Undyed)				
	5/8 Circle	27 mm	 Round Bodied	70 mm			SN 2473		

# TRUGLYDE








## Braided & Coated Polyglycolic Acid Suture

### Atraumatic Needled Sutures



NEEDLE				THREAD				
	Profile	Length	Description	Length	Size			
					3-0 (2)	2-0 (3)	0 (3.5)	1 (4)
	1/2 Circle	30 mm	 Taper Cutting	20 45 90 cm		SN 2511 AA SN 2511 A SN 2511	SN 2512 AA SN 2512 A SN 2512	SN 2513 AA SN 2513 A SN 2513
	1/2 Circle	30 mm	 Reverse Cutting	90 cm			SN 2535	
	1/2 Circle Double Armed	30 mm	 Rounded Bodied & Taper Cutting	140 cm			SN 2317 DA	
	1/2 Circle Heavy	30 mm	 Round Bodied	20 cm			SN 2344	SN 2343
	1/2 Circle	35 mm	 Taper Cutting	90 cm		SN 2762 F (Undyed)	SN 2763 F (Undyed)	
	1/2 Circle	36 mm	 Reverse Cutting OR Needle	90 cm			SN 2534	
	1/2 Circle	36 mm	 Taper Cutting	90 cm	SN 2516			
	1/2 Circle Double Armed	36 mm	 Reverse Cutting & Round Bodied	140 cm			SN 2777 F (Undyed)	
	1/2 Circle	40 mm	 Reverse Cutting OR Needle	90 cm				SN 2421
	1/2 Circle	40 mm	 Round Bodied	20 45 70 90 110 cm			SN 2346 AA SN 2346 A SN 2546 SN 2346 SN 2346 VS	
	1/2 Circle	40 mm	 Taper Cutting	90 110 cm		SN 2517	SN 2518 SN 2518 VS	
	1/2 Circle Double Armed	40 mm	 Rounded bodied & Taper Cutting	140 cm			SN 2346 DA	SN 2347 DA

# TRUGLYDE

## Braided & Coated Polyglycolic Acid Suture Atraumatic Needled Sutures

NEEDLE				THREAD				
	Profile	Length	Description	Length	Size			
					2-0 (3)	0 (3.5)	1 (4)	2(5)
	1/2 Circle	40 mm	Reverse Cutting	90 cm	SN 2382	SN 2422	SN 2350	
	1/2 Circle	40 mm HEAVY	Round Bodied	20 45 70 70 x 2 90 110			SN 2347 AA SN 2347 A SN 2359 SN 2359 DS SN 2347 SN 2347 VS	
	1/2 Circle	40 mm HEAVY	Taper Cutting	45 70 90 110		SN 2361	SN 2362 SN 2363 SN 2364 SN 2364 VS	
	1/2 Circle	45 mm	Round Bodied	90 cm				SN 2228
	3/8 Circle	45 mm	Round Bodied	70 cm			SN 2326	
	1/2 Circle	50 mm HEAVY	Reverse Cutting	90 cm			SN 2438	
	1/2 Circle	50 mm HEAVY	Round Bodied	100 cm			SN 2448	

## OB Kit (SN 2347 OB) contains one foil of each SN 2347DA & SN 2512A

	Profile	Length	Description	Length	Size			
					2-0 (3)	0 (3.5)	1 (4)	2(5)
	1/2 Circle Double Armed	40 mm	Round Bodied & Taper Cutting	140 cm			SN 2347 DA	
	1/2 Circle	30 mm	Taper Cutting	45 cm		SN 2512 A		