









TRUGUT

Sterilised Surgical Gut Suture: Plain

Atraumatic Needles and Sutures

























NEEDLE				THREAD			
	Profile	Length	Description	Length	Size		
				5-0(1.5)	4-0(2)	3-0(3)	2-0(3.5)
	3/8 Circle	12 mm	Reverse Cutting	76 cm	SN 4269		
	3/8 Circle	16 mm	Cutting	76 cm	SN 4270	SN 4271	
	3/8 Circle	22 mm	Cutting	76 cm	SN 5010		
	1/2 Circle	25 mm	Round Bodied	76 cm	SN 4238		
	3/8 Circle	26 mm	Cutting	76 cm	SN 5020		
	3/8 Circle	30 mm	Round Bodied	76 cm	SN 5033		
	Straight	50 mm	Round Bodied	76 cm	SN 5007	SN 5008	
	Straight	60 mm	Cutting	76 cm	SN 4288		

Standard and Short Length Sutures

Length	5-0 (1.5)	4-0 (2)	3-0 (3)	2-0 (3.5)	0 (4)	1 (5)	2 (6)	3 (7)	4 (8)
1 x 152 cm	S 2000	S 2001	S 2002	S 2003	S 2004	S 2005	S 2006	S 2007	S 2008
2 x 76 cm	S 2240	S 2241	S 2242	S 2243	S 2244	S 2245	S 2246	S 2247	S 2248

TRUGUT












Sterilised Surgical Gut Suture: Chronic Atraumatic Needles and Sutures

NEEDLE				THREAD					
	Profile	Length	Description	Length	5-0(1.5)	4-0(2)	Size 3-0(3)	2-0(3.5)	0 (4)
	3/8 Circle	12 mm	 Reverse Cutting	38 76 cm			SN 4266		
	1/2 Circle	12 mm	 Reverse Cutting	76 cm				SN 4272	
	3/8 Circle	16 mm	 Round Bodied	76 cm	SN 4281	SN 5048			
	3/8 Circle	16 mm	 Cutting	76 cm			SN 4280	SN 4274	
	5/8 Circle	16 mm	 Cutting	38 cm				SN 4265	
	1/2 Circle	20 mm	 Round Bodied	76 cm				SN 4237	
	3/8 Circle	22 mm	 Cutting	76 cm				SN 4201	
	1/2 Circle	26 mm	 Round Bodied	76 cm			SN 4125	SN 4123	SN 4124
	3/8 Circle	26 mm	 Reverse Cutting	76 cm				SN 4828	
	3/8 Circle	30 mm	 Round Bodied	76 cm				SN 4216	SN 4217
	1/2 Circle	30 mm	 Round Bodied	76 152 cm				SN 4241	SN 4242 SN 2214
	1/2 Circle	30 mm HEAVY	 Round Bodied	76 cm					SN 4254
	1/2 Circle	30 mm HEAVY	 Cutting	76 cm				SN 4255	
	3/8 Circle	32 mm	 Round Bodied	76 cm				SN 4113	SN 4114

TRUGUT

Sterilised Surgical Gut Suture: Chromic











Atraumatic Needles and Sutures

NEEDLE				THREAD				
Profile	Length	Description	Length	Size				
				3-0 (3)	2-0 (3.5)	0 (4)	1 (5)	2 (6)
	5/8 Circle 35 mm	Round Bodied	76 cm	SN 4264				
	1/2 Circle 36 mm Double Armed	1/2 Circle Round Bodied 1/2 Circle Reverse Cutting	152 cm	SN 4277				
	1/2 Circle 38 mm	Round Bodied	76 cm	SN 4295	SN 4296	SN 4298	SN 4299	
	Straight 40 mm	Round Bodied	76 cm	SN 4200				
	3/8 Circle 40 mm	Round Bodied	76 cm	SN 4220	SN 4221			
	1/2 Circle 40 mm	Round Bodied	76 cm 75 cm	SN 4294	SN 4245	SN 4246		SN 4260
	1/2 Circle 40 mm HEAVY	Round Bodied	76 cm 152 cm			SN 4258	SN 4259 SN 4221	
	1/2 Circle 40 mm HFAVY	Trocar Point	100 cm 76 cm				SN 4239	SN 4244 SN 4239A
	1/2 Circle 40 mm HEAVY	Taper Cutting	100 cm 76 cm			SN 4279 SN 4257	SN 4297 SN 4297A	
	Straight 44 mm	Round Bodied	76 cm	SN 4102				
	Straight 45 mm	Round Bodied	76 cm	SN 4204	SN 4205			

TRUGUT

Sterilised Surgical Gut Suture: Chromic









Atraumatic Needles and Sutures

NEEDLE			THREAD				
Profile	Length	Description	Length	Size			
			3-0 (3)	2-0 (3.5)	0 (4)	1 (5)	2 (6)
 3/8 Circle	45 mm	Round Bodied	76 cm	SN 4223	SN 4224	SN 4225	SN 4226
 3/8 Circle Double Armed	45 mm	Round Bodied	76 cm				SN 4826
 1/2 Circle	45 mm	Round Bodied	76 cm		SN 4249	SN 4250	
 1/2 Circle Double Armed	45 mm	1/2 Circle Round Bodied and 1/2 Circle Reverse Cutting	152 cm			SN 4276	SN 4275
 1/2 Circle	45 mm HEAVY	Round Bodied	100 cm			SN 4229	SN 4227 SN 4228
 Straight	50 mm	Round Bodied	76 cm	SN 4207	SN 4208	SN 4209	SN 4210
 1/2 Circle	50 mm HEAVY	Round Bodied	100 cm				SN 4248
 1/2 Circle	50 mm	Reverse Cutting	76 cm				SN 4273
 Straight	60 mm	Round Bodied	76 cm		SN 4212	SN 4213	
 3/8 Circle	60 mm	Cutting	76 cm				SN 4262

TRUGUT

Sterilised Surgical Gut Suture: Chromic

Atraumatic Needles and Sutures

NEEDLE			THREAD			
Profile	Length	Description	Length	Size		
 5/8 Circle	60 mm	 Round Bodied	76 cm	2-0 (3.5)	0 (4)	1 (5) 2 (6)
 3/8 Circle BLUNT	63 mm	 Round Bodied	76 cm	SN 4267		
 3/8 Circle	65 mm	 Round Bodied	76 cm	SN 4234	SN 4233	
 1/2 Circle	65 mm	 Round Bodied	100 cm	SN 4278		

Standard and Short Length Sutures

LENGTH	SIZE							
	5-0 (1.5)	4-0 (2)	3-0 (3)	2-0 (3.5)	0 (4)	1 (5)	2 (6)	3(7)
1 x 152 cm	S 2210	S 2211	S 2212	S 2213	S 2214	S 2215	S 2216	S 2217
2 x 76 cm	S 2250	S 2251	S 2252	S 2253	S 2254	S 2255	S 2256	

TRUGUT

Technical Details

- **TRUGUT Catgut Sutures**
Absorbable Surgical Suture, USP
Catgut Plain & Chromic
- **Raw Material Used**
100% submucosal layer of sheep
intestine.
- **Colour**
Brown(Chromic), Yellow (Plain)
- **Dye Used (Chromic)**
Pyrogallol (US FDA approved colour)
- **Sterilization**
Ethylene Oxide
- **Shelf Life**
5 years
- **Range**
USP Sizes available - 5-0 to 2.
Standard sutures - 150cm, 2x76cm
With needle Suture lengths from
38cm to 150cm & Needle length from
12mm to 65mm
- **Packing**
Wet (Packed in Iso propyl alcohol-IPA)
Primary packaging
Aluminium foil
Secondary packaging
Medical grade paper & Polyester
laminated poly film.
Pack of 12pcs
400gsm Super Chromo board.
Master shipper
100 boxes of 12s (5 ply of 20 Kgs/cm²
burst strength - 54 x 33 x 33 cm)
weighing approximately 10 kgs.

Features

- Absorption within 60-90
days (Chromic).
Absorption within 45-60
days (Plain)
- Highest level of collagen
guarantees predictable
absorption.
- Accurate polishing ensures
smooth tissue passage.
- Packed in IPA to reduce
memory & increase
pliability.
- Ribbons individually
chromicised to assure
uniform absorption &
wound support.
- Rigidly controlled spinning
& drying to ensure higher
tensile strength



SUTURES INDIA PRIVATE LIMITED

472 D, 13th Cross, 4th Phase, Peenya Industrial Area, Bangalore-560 058. INDIA
Tel : +91 80 4186 8000 (30 Lines) ; Fax : +91 80 4117 1056, Email : sales@suturesin.com
www.suturesin.com