

Monofilament Polypropylene Suture

TRULENE

Atraumatic Needles and Sutures



USP 8-0	Needle Profile	Needle Description	Needle Length	Suture Length	Suture Code
	3/8 Circle	▼ Taper Point Double Armed	6 mm	60 cm	SN 2086
	3/8 Circle	▼ Taper Point Double Armed	6.5 mm	60 cm	SN 8732
	3/8 Circle	▼ Taper Point Double Armed	8 mm	45 cm	SN 2087
	3/8 Circle	▼ Taper Point Double Armed	8 mm	60 cm	SN 2088
	3/8 Circle	▼ Taper Point Double Armed	8 mm	75 cm	SN 2088 L
	3/8 Circle	▼ Antiglare Taper Point Double Armed	8 mm	60 cm	SN 8614 B
	3/8 Circle	▼ Taper Point Double Armed	9 mm	60 cm	SN 8741
USP 7-0					
	3/8 Circle	▼ Taper Point Double Armed	6 mm	60 cm	SN 2076
	3/8 Circle	▼ Taper Point Double Armed	6.5 mm	60 cm	SN 8731
	3/8 Circle	▼ Taper Point Double Armed	7.6 mm	60 cm	SN 8735
	3/8 Circle	▼ Taper Point Double Armed	8 mm	45 cm	SN 2078 S
	3/8 Circle	▼ Taper Point Double Armed	8 mm	60 cm	SN 2078
	3/8 Circle	▼ Antiglare Taper Point Double Armed	8 mm	60 cm	SN 8304 B
	3/8 Circle	▼ Taper Point Double Armed	9 mm	60 cm	SN 2079
	3/8 Circle	▼ Antiglare Taper Point Double Armed (BP)	9 mm	60 cm	SN 8121 B
	3/8 Circle	▼ Taper Point Double Armed	9 mm	75 cm	SN 2079 L
	3/8 Circle	▼ Taper Point Double Armed	9.7 mm	60 cm	SN 8702
	3/8 Circle	▼ Taper Point Double Armed	10 mm	60 cm	SN 2070
USP 6-0					
	3/8 Circle	▼ Antiglare Taper Point Double Armed (BP)	8 mm	60 cm	SN 8305 B
	3/8 Circle	▼ Taper Point Double Armed	9 mm	45 cm	SN 8806
	3/8 Circle	▼ Taper Point Double Armed	9 mm	60 cm	SN 8805
	3/8 Circle	▼ Taper Point Double Armed	9.7 mm	60 cm	SN 8805 M
	1/2 Circle	▼ Taper Point Double Armed	10 mm	60 cm	SN 8207
	1/2 Circle	▼ Taper Point Double Armed (Pledget 3mm x 3mm x 1.5mm)	10 mm	60 cm	SN 8207 PL 3
	3/8 Circle	▼ Taper Point Double Armed	10 mm	60 cm	SN 8804
	3/8 Circle	▼ Taper Point Double Armed	10 mm	75 cm	SN 8804 L
	1/2 Circle	▼ Taper Point Double Armed	10 mm	90 cm	SN 8611
	3/8 Circle	▼ Reverse Cutting Double Armed	12 mm	90 cm	SN 8612



SUTURES INDIA PRIVATE LIMITED

472 D, 13th Cross, 4th Phase, Peenya Industrial Area, Bangalore-560 058. INDIA

Tel : +91 80 4186 8000 (30 Lines); Fax : +91 80 4117 1056,

Email : sales@suturesin.com; www.suturesin.com

USP 6-0	Needle Profile	Needle Description	Needle Length	Suture Length	Suture Code
	3/8 Circle	▼ Taper Point Double Armed	12.3 mm	60 cm	SN 8726 M
	3/8 Circle	▼ Taper Point Double Armed	13 mm	45 cm	SN 8726 S
	3/8 Circle	▼ Taper Point Double Armed	13 mm	60 cm	SN 8726
	1/2 Circle	▼ Taper Point Double Armed	13 mm	60 cm	SN 8727
	3/8 Circle	▼ Antiglare Taper Point Double Armed	13 mm	60 cm	SN 8307 B
	3/8 Circle	▼ Antiglare Taper Point Double Armed	13 mm	75 cm	SN 8725 B
	1/2 Circle	▼ Taper Point Double Armed	13 mm	90 cm	SN 8613
USP 5-0					
	1/2 Circle	▼ Taper Point Double Armed	10 mm	90 cm	SN 8610
	½ Circle	▼ Antiglare Taper Point Double Armed	11 mm	75 cm	CJ 8802
	1/2 Circle	▼ Taper Point Double Armed	13 mm	45 cm	SN 8710 S
	1/2 Circle	▼ Taper Point Double Armed	13 mm	60 cm	SN 8205
	3/8 Circle	▼ Antiglare Taper Point Double Armed	13 mm	60 cm	SN 8321 B
	1/2 Circle	▼ Taper Point Double Armed	13 mm	75 cm	SN 8710
	1/2 Circle	▼ Taper Point Double Armed (Pledget 3mm x 3mm x 1.5mm)	13 mm	75 cm	SN 8710 PL 3
	1/2 Circle	▼ Taper Point Double Armed	16 mm	90 cm	SN 8556
	1/2 Circle	▼ Taper Point Double Armed (Pledget 3mm x 3mm x 1.5mm)	16 mm	90 cm	SN 8556 PL 3
USP 4-0					
	1/2 Circle	▼ Taper Point Double Armed	13 mm	60 cm	SN 8204
	1/2 Circle	▼ Taper Point Double Armed	13 mm	60 cm	SN 8204 PL 3
	1/2 Circle	▼ Taper Point Double Armed	16 mm	90 cm	SN 8557
	1/2 Circle	▼ Taper Cut Double Armed	17.5 mm	90 cm	SN 8935
	1/2 Circle	▼ Taper Cut Double Armed (Pledget 3mm x 3mm x 1.5mm)	17.5 mm	90 cm	SN 8935 PL 3
	1/2 Circle	▼ Taper Point Double Armed	26 mm	90 cm	SN 8521
USP 3-0					
	1/2 Circle	▼ Taper Cut Double Armed	17.5 mm	90 cm	SN 8936
	1/2 Circle	▼ Taper Point	26 mm	75 cm	SN 8832
	1/2 Circle	▼ Taper Cut Double Armed	26 mm	90 cm	SN 8976
	1/2 Circle	▼ Taper Point Double Armed	26 mm	90 cm	SN 8522
USP 2-0					
	1/2 Circle	▼ Taper Cut Double Armed	17.5 mm	90 cm	SN 8937
	1/2 Circle	▼ Taper Cut Double Armed	26 mm	90 cm	SN 8977
	1/2 Circle	▼ Taper Cut Double Armed	26 mm	150 cm	SN 8977 L
	1/2 Circle	▼ Taper Point Double Armed	30 mm	90 cm	SN 8526